



GIO.ZIRC

OUR GIO.ZIRC ZIRCONIA CROWN IS AN EXCEPTIONAL RESTORATION WITH STUNNING AESTHETICS AND STRENGTH LAYERED WITH NORITAKE PORCELAIN

 WHITE PEAKS
DENTAL SOLUTIONS

Noritake

GIOGIO
dental laboratory

we're all smiles.®

Produced utilising our in-house, Industry leading CAD-CAM facilities with premium Ultra-translucent zirconium oxide by White Peaks.

This material is fully biocompatible and extremely robust with over 65% higher strength than lithium disilicate glass ceramics.


We hand layer chip resistant Noritake porcelain which is specifically produced for Zirconia substructures. The results are exceptionally beautiful and undetectable.

Our Manufacturing process ensures perfect fit and predictable results. It can be used to make single copings and large bridges as well as various indications in Implantology.


IMPRESSION TECHNIQUE	Digital or Traditional
MATERIAL	Zirconium Dioxide layered with Porcelain
RECOMMENDED USE	Single Units Crowns, Bridges
AESTHETICS	*****
CEMENTATION	Regular
MANUFACTURING METHOD	CAD-CAM milled coping with hand layered porcelain
STRENGTH	*****
WARRANTY	5 years
TURNAROUND TIME	5 working days
MARGINS	Feather edge, shoulder or chamfer
HIGHLIGHTED BENEFITS	Metal free biocompatibility, High strength and aesthetics



PREPARATION GUIDE



Recommended Preparation For Anterior Teeth



Recommended Preparation For Posterior Teeth