



IPS EMPRESS®

OVER 20 YEARS OF TRIED AND TESTED CLINICAL EXPERIENCE WITH LEUCITE GLASS CERAMIC THAT HAS SET DECISIVE STANDARDS IN TRANSLUCENCY AND LIFELIKENESS. FORTY MILLION RESTORATIONS PRODUCED AND COUNTING.

IPS **Empress**®

GIORGIO
dental laboratory

we're all smiles.®

The IPS Empress® glass-ceramic is an extremely homogeneous, leucite-based material that scatters light naturally and provides a balanced chameleon effect for optimum shade matching.


Originally, like e.max®, Empress® could only be processed by waxing and pressing however since the introduction of IPS Empress® CAD, the popular all-ceramic material may now additionally be processed in a modern fashion using the CAD/CAM technology.

Giorgio Dental has offered this restoration since the early 90s and we continue to use it for the many different indications available because of its accurate predictability, it's beautiful aesthetics and its relative ease of use for Dentists and Technicians alike.


IMPRESSION TECHNIQUE	Digital or Traditional
MATERIAL	Leucite Glass Ceramic
RECOMMENDED USE	Single tooth, inlays/onlays veneers, anteriors and posteriors
AESTHETICS	*****
CEMENTATION	Regular or adhesive bonding
MANUFACTURING METHOD	CAD-CAM milled or waxed and pressed by hand
STRENGTH	****
WARRANTY	5 Years on crowns
TURNAROUND TIME	5 working days
MARGINS	Shoulder or chamfer
HIGHLIGHTED BENEFITS	Highest Aesthetics, strength, metal-free

IPS Empress®

PREPARATION GUIDE



Recommended Preparation For Anterior Teeth



Recommended Preparation For Posterior Teeth