



PFM

(PORCELAIN FUSED TO METAL)

DESPITE THE QUICKENING DEVELOPMENT AND INTRODUCTION OF NEW MATERIALS IN THE INDUSTRY, THE DEPENDABLE AND EVER-POPULAR PFM RESTORATION IS A QUALITY OPTION FOR ANY DENTIST SEEKING BEAUTY AND BRAWN.

 WHITE PEAKS
DENTAL SOLUTIONS

GIORGIO
dental laboratory

we're all smiles.®

The PFM has been the gold standard and star of Dental Restorations since its invention in the late 1950s by Abraham Weinstein.

Whilst the method of manufacture and the quality of materials have changed significantly since those days, the attention-to-detail and standard of workmanship applied to our traditional PFM remains the same.


The humble and widely used PFM can be constructed with a range of materials including non-precious nickel free chrome cobalt and precious gold alloy.

It is still a popular choice amongst aesthetically and restoratively minded Dentists everywhere due to its combination of marginal integrity, strength, biocompatibility and aesthetics.


IMPRESSION TECHNIQUE	Digital or Traditional
MATERIAL	Non-precious or Precious Gold Alloy
RECOMMENDED USE	Full Mouth, single units, Implants
AESTHETICS	★★★★
CEMENTATION	Regular
MANUFACTURING METHOD	CAD-CAM milled coping with hand layered porcelain
STRENGTH	★★★★
WARRANTY	5 years
TURNAROUND TIME	5 working days
MARGINS	Feather edge, shoulder or chamfer
HIGHLIGHTED BENEFITS	Traditional restoration, alloy choice, layered build



PREPARATION GUIDE



Recommended Preparation For Anterior Teeth



Recommended Preparation For Posterior Teeth