



VITA SUPRINITY®

VITA DELIVERS STUNNING AESTHETICS BASED ON INTEGRATED TRANSLUCENCY, OPALESCENCE AND FLUORESCENCE. THE INNOVATIVE GLASS CERAMIC IS ENRICHED WITH ZIRCONIA RESULTING IN A HIGH-STRENGTH, ZIRCONIA-REINFORCED LITHIUM SILICATE CERAMIC (ZLS).

VITA

GIORGIO
dental laboratory

we're all smiles.®

Whilst we are happy to accept traditional impressions, VITA Suprinity® is now suitable for use within our Digital Workflow.

Users of Suprinity® indicate high reliability thanks to zirconia-reinforced glass ceramic matrix, optimized edge stability and mechanical load capacity leading to long term clinical success.


It achieves brilliant aesthetics based on integrated translucency, opalescence and fluorescence.

This incredibly versatile material can be applied to many different indications including Crowns, Veneers, inlays/onlays and implants.


IMPRESSION TECHNIQUE	Digital or Traditional
MATERIAL	Zirconia-reinforced lithium silicate ceramic (ZLS)
RECOMMENDED USE	Full mouth, veneers, onlays/inlays, implants
AESTHETICS	*****
CEMENTATION	Regular or adhesive bonding
MANUFACTURING METHOD	CAD-CAM milled coping with hand layered porcelain
STRENGTH	*****
WARRANTY	5 years on crowns
TURNAROUND TIME	5 working days
MARGINS	Feather edge, shoulder or chamfer
HIGHLIGHTED BENEFITS	Highest Aesthetics, strength, metal-free

VITA

PREPARATION GUIDE



Recommended Preparation For Anterior Teeth



Recommended Preparation For Posterior Teeth